


DIAGNOSTICS INSTRUMENTS

MCINTOSH-LARYNGOSCOPE

MILLER-LARYNGOSCOPE

LARYNGOSCOPE INSTRUMENTS



MANUFACTURER OF HIGH-QUALITY DIAGNOSTICS INSTRUMENTS

Mc – Intosh Laryngoscopes Standard Illumination Adult Size



SP-8904

Handle Standard Illumination
 Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 0, 1, 2 and Adult Handle	▪ SP-8905
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ SP-8906
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ SP-8907
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ SP-8907

Sets available in hard case or soft pouch



SP-8900

Mc - Intosh Blade Standard Illumination
 Size: 0



SP-8901

Mc - Intosh Blade Standard Illumination
 Size: 1



SP-8902

Mc - Intosh Blade Standard Illumination
 Size: 2



SP-8903

Mc - Intosh Blade Standard Illumination
 Size: 3



SP-8908

Mc - Intosh Blade Standard Illumination
 Size: 4

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7276 Standard



Mc – Intosh Laryngoscopes

Fiber Optic Illumination

Adult Size



SP-8913

Handle Fiber Optic Illumination
Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 0, 1, 2 and Adult Handle	■ 21-200-30
■ Set of 3 Blades Size 1, 2, 3 and Adult Handle	■ 21-200-31
■ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	■ 21-200-40
■ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	■ 21-200-41

Sets available in hard case or soft pouch



SP-8909

Mc - Intosh Blade Fiber Optic Illumination
Size: 0



SP-8910

Mc - Intosh Blade Fiber Optic Illumination
Size: 1



SP-8911

Mc - Intosh Blade Fiber Optic Illumination
Size: 2



SP-8912

Mc - Intosh Blade Fiber Optic Illumination
Size: 3



SP-8914

Mc - Intosh Blade Fiber Optic Illumination
Size: 4

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel



Miller Laryngoscopes Standard Illumination Child / Infant Size



SP-8919

Handle Standard Illumination
Size: Child / Infant

Handle options: Adult & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 00, 0, 1 and Child Handle	■ 22-200-30
■ Set of 3 Blades Size 0, 1, 2 and Child Handle	■ 22-200-31
■ Set of 4 Blades Size 00, 0, 1, 2 and Child Handle	■ 22-200-40
■ Set of 4 Blades Size 0, 1, 2, 3 and Child Handle	■ 22-200-41

Sets available in hard case or soft pouch



SP-8915

Miller Blade Standard Illumination
Size: 00



SP-8916

Miller Blade Standard Illumination
Size: 0



SP-8917

Miller Blade Standard Illumination
Size: 1



SP-8918

Miller Blade Standard Illumination
Size: 2



SP-8920

Miller Blade Standard Illumination
Size: 3

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination



Miller Laryngoscopes Fiber Optic Illumination Adult Size



SP-8925

Handle Fiber Optic Illumination
Size: Adult
Handle options: Child & Stubby

Standard Sets

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	23-200-30
Set of 3 Blades Size 1, 2, 3 and Adult Handle	23-200-31
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	23-200-40
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	23-200-41

Sets available in hard case or soft pouch



SP-8921

Miller Blade Fiber Optic Illumination
Size: 0



SP-8922

Miller Blade Fiber Optic Illumination
Size: 1



SP-8923

Miller Blade Fiber Optic Illumination
Size: 2



SP-8924

Miller Blade Fiber Optic Illumination
Size: 3



SP-8926

Miller Blade Fiber Optic Illumination
Size: 4

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable



Mc – Intosh Integrated Laryngoscopes Fiber Optic Illumination Adult Size



SP-8933

Handle Fiber Optic Illumination
Size: Adult

Handle options: Child & Stubby

Standard Sets

Contents	Code
■ Set of 3 Blades Size 00, 0, 1 and Adult Handle	■ 24-200-30
■ Set of 3 Blades Size 1, 2, 3 and Adult Handle	■ 24-200-31
■ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	■ 24-200-40
■ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	■ 24-200-41

Sets available in hard case or soft pouch



SP-8927

Mc - Intosh Blade Fiber Optic Illumination
Size: 00



SP-8928

Mc - Intosh Blade Fiber Optic Illumination
Size: 0



SP-8929

Mc - Intosh Blade Fiber Optic Illumination
Size: 1



SP-8930

Mc - Intosh Blade Fiber Optic Illumination
Size: 2



SP-8931

Mc - Intosh Blade Fiber Optic Illumination
Size: 3



SP-8932

Mc - Intosh Blade Fiber Optic Illumination
Size: 4

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable



McCoy Flexi - Tip Laryngoscopes

Fiber Optic Illumination

Adult Size



SP-8937

Handle Fiber Optic Illumination
Size: Adult



SP-8934

McCoy Blade Fiber Optic Illumination
Size: 1



SP-8935

McCoy Blade Fiber Optic Illumination
Size: 2



SP-8936

McCoy Blade Fiber Optic Illumination
Size: 3



Standard Sets

Contents	Code
Set of 2 Blades Size 2, 3 and Adult Handle	25-200-30
Set of 2 Blades Size 3, 4 and Adult Handle	25-200-31
Set of 3 Blades Size 2, 3, 4 and Adult Handle	25-200-40

Sets available in hard case or soft pouch

Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel



McCoy Flexi - Tip Integrated Laryngoscopes Fiber Optic Illumination Adult Size



SP-8941

Handle Fiber Optic Illumination
Size: Adult



SP-8938

McCoy Blade Fiber Optic Illumination
Size: 1



SP-8939

McCoy Blade Fiber Optic Illumination
Size: 2



SP-8940

McCoy Blade Fiber Optic Illumination
Size: 3

Standard Sets

Contents	Code
■ Set of 2 Blades Size 2, 3 and Adult Handle	■ 26-200-30
■ Set of 2 Blades Size 3, 4 and Adult Handle	■ 26-200-31
■ Set of 3 Blades Size 2, 3, 4 and Adult Handle	■ 26-200-40

Sets available in hard case or soft pouch



Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable





SP-8942

Mini Otoscope Fiber Optic
Complete Set of Disposable Cannulas



Color Options

Specification:

- Manufactured From Plastic Material
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V





Fiber Optic Diagnostic Set

SP-8943



Crafted with Care
Small Diagnostic Set

SP-894

Universal Diagnostic Set

SP-8945





SP-8946

Brass Otoscope Fiber Optic
Brass Otoscope Conventional



SP-8947

ENT Set Complete



SP-8948

Brass Ophthalmoscope

Specification:

- Manufactured From Brass Chrome plated
- 3x Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries





Baby Prestige stethoscope

SP-8949



Child Prestige stethoscope

SP-8950



Kirchner Colorscope stethoscope

SP-8951



Standard stethoscope

SP-8952



Dual Head stethoscope

SP-8953



Standard Accessories

Nasal Speculum



SP-8954



SP-8955



SP-8956

Bulbs



SP-8957



SP-8958



SP-8959



SP-8960

Ear Speculas



SP-8961



SP-8962



SP-8963



SP-8964

Ear Speculas & Testing Tube



SP-8965



SP-8966



SP-8967



SP-8968



SP-8969



SP-8970



SP-8971



Company Name: SurgiPortal

Phone: +92 339 0059380

Email: info@surgiportal.com

Location: Ikram Plaza, Near KFC, Kashmir Road, Sialkot, Pakistan