

EYE SURGERY INSTRUMENTS

BASIC EYE SURGERY ITEMS

ALL EYE SURGERY SETS

EYE SURGERY FORCEPS



MANUFACTURER OF EYE SURGERY INSTRUMENTS

INDEX

FORCEPS

Tying	3
Toothed	4
Notched	5
Chalazion	5
Capsulorhexis	6
Lens	7
Hemostatic	8
Prechopper	8

NEEDLE HOLDERS 9

SCISSORS 10

HOOKS 13

CHOPPERS 14

FIXATION RING 14

MANIPULATORS 15

SPATULA 15

SPECULUMS 16

CALIPERS 17

PUNCHES 18

LASIK 18

VITREORETINAL 19

IRRIGATION / ASPIRATION 25


SETS 26

STERILIZATION TRAYS 27

2. TOOTHED FORCEPS

Straight Toothed Forceps

1 x 2 teeth, with 5.0mm tying platforms




SP-8233
SP-8233
SP-8233
SP-8233



85mm long | 0.10mm teeth
85mm long | 0.12mm teeth
115mm long | 0.10mm teeth
115mm long | 0.12mm teeth

Colibri Conjunctiva Toothed Forceps

1 x 2 teeth, with 5.0mm tying platforms




SP-8233



85mm long | 0.12mm teeth

Bonn Toothed Forceps

1 x 2 teeth



SP-8233



85mm long | 0.12mm teeth

Colibri Toothed Forceps

1 x 2 teeth



SP-8233
SP-8233



85mm long | 0.12mm teeth with tying platform
115mm long | 0.12mm teeth with tying platform

Harms-Colibri Toothed Forceps

1 x 2 teeth | 75mm long



SP-8233
SP-8233
SP-8233




0.12mm teeth with tying platform
0.30mm teeth with tying platform
0.50mm teeth with tying platform

FORCEPS

1. TYING FORCEPS

Straight Tying Forceps

Straight shafts with 4.50mm tying platforms



SP-8233

SP-8233

SP-8233

SP-8233



85mm long | Standard tips

85mm long | Ultrafine tips

115mm long | Standard tips

115mm long | Ultrafine tips

Curved Tying Forceps

Curved shafts with 6mm tying platforms | ideal for 10/0 or 11/0 suture



SP-8233

SP-8233




115mm long

85mm long

McPherson Tying Forceps

Angled 4.50mm tying platforms



SP-8233

SP-8233

SP-8233

SP-8233

SP-8233

SP-8233



85mm long | Standard tips | 12mm angled shafts

85mm long | Ultrafine tips | 7mm angled shafts

85mm long | Ultrafine tips | 10mm angled shafts

85mm long | Ultrafine tips | 12mm angled shafts

115mm long | Standard tips

115mm long | Ultrafine tips

3. NOTCHED FORCEPS

Straight Notched Forceps

0.25mm notched | with tying platforms | 85mm long



SP-8233
SP-8233



Standard Tips
Ultrafine Tips

Angled Notched Forceps

0.25mm notched | with tying platforms | 85mm long



SP-8233
SP-8233



Standard Tips
Ultrafine Tips

Pierse Corneal Forceps

0.25 mm notched | very delicate | 115 mm long



SP-8233
SP-8233



Straight shafts
Curved shafts

4. CHALAZION FORCEPS

Hunt Chalazion Forceps

98mm long



SP-8233
SP-8233
SP-8233
SP-8233



Oval 8x11mm
Round 8mm
Oval 12x14 mm
Round 15mm

5. CAPSULORHEXIS FORCEPS

Utrata-Style Capsulorhexis Forceps

Fine angled sharp tips, 12mm 45°angled shafts



Utrata-Style Capsulorhexis Forceps, *curved*

Fine angled sharp tips, 12mm 45°angled and curved shafts



Utrata-Style Capsulorhexis Forceps

45° angled tips with slim 10mm shafts | Round handle for fingertip rotational control



105 mm long | With guide pin | Iris stop |

Corydon-Style Capsulorhexis Forceps, cystotome teeth

Fine 45°angled sharp tips, with curved 12mm shafts



85mm long
115mm long

Corydon-Style Capsulorhexis Forceps, cystotome teeth

Fine 45° angled tips with curved 12 mm shafts | Round handle for fingertip rotational control




105 mm long | With guide pin | Iris stop

7. HEMOSTATIC FORCEPS

Hartmann Mosquito Forceps

Ring handle with ratchet lock, 95mm long



SP-8233

SP-8233



Straight 20mm serrated jaws
Curved 20mm long serrated jaws

8. PRECHOPPER

Akahoshi PreChopper

Conventional smooth sharp pointed tips | 12mm angled curved shafts | cross-action



SP-8233



3 mm maximum tip opening | 120mm long

Lehner-Utrata Style Capsulorhexis Forceps

Fine angled sharp tips with curved 12mm shafts for capsulorhexis through tunnel incisions



SP-8233



105 mm long | With guide pin

Capsulorhexis Forceps

45° angled and Iris stop



SP-8233

SP-8233



12mm straight shafts
12mm curved shafts

6. LENS FORCEPS

Kelman-McPherson Lens Forceps

7.5mm long smooth jaws | very delicate



SP-8233



108mm long

Shepard Lens Forceps

Gently curved shafts with flat Tungsten Carbide coated jaws



SP-8233



125mm long

NEEDLE HOLDERS

Barraquer Needle Holder - Curved

Delicate, 9mm gently curved jaws, sharp tips, 115mm long



SP-8233

SP-8233



With Lock

Without Lock

Troutman Needle Holder

Curved, 6mm delicate jaws, 110mm long



SP-8233

SP-8233



With Lock

Without Lock

SCISSORS

Conjunctiva Scissors

Gently curved blades, 12mm tip to pivot length, blunt tips



SP-8233

115mm long

Eye Scissors

Straight blade, Pointed tips, 32mm from pivot to tip



SP-8233

110mm long

Castroviejo Corneal Section Scissors

Blunt tips, 0.5mm longer lower blades than upper ones, 9mm from pivot to tip, 105mm long



SP-8233

Right
Left

SP-8233

Iris Scissors

Ring handle, Pointed tips



SP-8233

SP-8233

SP-8233

SP-8233

SP-8233

SP-8233

SP-8233

SP-8233

Straight tips, 18mm from mid-screw to tip, 98mm long
Curved tips, 18mm from mid-screw to tip, 98mm long
Straight Tips, 25mm from mid-screw to tip, 105mm long
Curved Tips, 25mm from mid-screw to tip, 105mm long
Straight Tips, 28mm from mid-screw to tip, 100mm long
Curved Tips, 28mm from mid-screw to tip, 100mm long
Sharp tips, 11mm blade length, Straight blades, 105mm long

Stitch Scissors

Ring handle, 24mm blade, 28mm from pivot to tip, 95mm long



SP-8233
SP-8233
SP-8233
SP-8233
SP-8233
SP-8233

Straight, pointed tips
Straight, pointed left tip, blunt right tip
Straight, blunt tips
Curved, pointed tips
Curved, pointed left tip, blunt right tip
Curved, blunt tips

Stitch Scissors

Ring handle, 26mm blade, 30mm from pivot to tip, 106mm long



SP-8233
SP-8233
SP-8233
SP-8233
SP-8233
SP-8233

Straight, pointed tips
Straight, pointed left tip, blunt right tip
Straight, blunt tips
Curved, pointed tips
Curved, pointed left tip, blunt right tip
Curved, blunt tips

Stitch Scissors

Ring handle, 15mm blade, 20mm from pivot to tip, 88mm long



SP-8233
SP-8233
SP-8233
SP-8233

Straight, pointed tips
Straight, blunt tips
Curved, pointed tips
Curved, blunt tips

Westcott Tenotomy Scissors



SP-8233
 SP-8233
 SP-8233
 SP-8233

19mm from pivot to tip, curved blunt tip, 110mm long
 21mm from pivot to tip, curved blunt tip, 115mm long
 23mm from pivot to tip, curved blunt tip, 140mm long
 23mm from pivot to tip, curved sharp tip, 140mm long

Vannas Capsulotomy Scissors

Sharp pointed tips, curved blades, 95mm long



SP-8233
 SP-8233

10.5mm blade
 12.5mm blade

Vannas Scissors



SP-8233
 SP-8233
 SP-8233
 SP-8233

Sharp pointed, straight 7.5mm blades, 9.5mm from pivot to tip, 82mm
 Sharp pointed, curved 7.5mm blades, 9.5mm from pivot to tip, 82mm
 Sharp pointed, side angled 7.0mm blades, 9.0mm from pivot to tip, 82mm
 Extra thin 7.0mm long blades, sharp pointed, 9.0mm from pivot to tip, 105mm

Osher IOL Cutter



SP-8233
 SP-8233

Side-angled
 Straight

CHOPPERS

Phaco Chopper



SP-8233
SP-8233



Straight
Angled

Phaco Chopper



SP-8233
SP-8233



Right hand
Left hand

Phaco Chopper

Double-ended



SP-8233



Paddle manipulator / 60° phaco chopper



Phaco Chopper

Double-ended for vertical and horizontal phaco chopping



SP-8233



Lieberman Micro Finger / Phaco Chopper



FIXATION RING

Fine/Thornton Fixation Ring

Swivel head




SP-8233
SP-8233

13mm/ 9mm cut out
16mm/14mm cut out

HOOKS

Kuglen Iris Hook and Lens Manipulator



SP-8233
SP-8233



Angled
Straight

Sinskey Lens Manipulating Hook



SP-8233



0.15 mm diameter, blunt tip, angled to right

Jameson Muscle Hook



SP-8233
SP-8233



1.5mm bulbous tips
2.0mm bulbous tips with 6mm flat portion

Graefe Strabismus Hook



SP-8233
SP-8233
SP-8233
SP-8233




Size 8mm x 1.5mm, 135mm long
Size 9mm x 1.5mm, 135mm long
Size 10mm x 1.5mm, 135mm long
Size 10mm x 2.0mm, 135mm long

SPECULUMS

Lindstrom-Chu Aspirating Speculum

Prevents excess pooling of BSS from eyelid margins. Unique reversible speculum design to rest the device nasally or temporally.



SP-8233

65mm long

Lieberman Speculum, K-wire


Adjustable wire speculum



SP-8233

45° angled shaft. 14.5mm open blades. 65mm long

Kratz-Barraquer Wire Speculum



SP-8233

13.5mm open blades, 52mm long

Barraquer Solid Speculum



SP-8233

14.5mm long blades, 52mm long

MANIPULATORS

Mushroom Nucleus Manipulator



SP-8233



Drysdale Nucleus Manipulator



SP-8233



Paddle-shaped tip

Castroviejo Cyclodialysis/Synechia Spatula

Double-ended

SP-8233

SP-8233

0.25mm wide
0.50mm wide



Membrane Spatula

SP-8233

SP-8233

SP-8233


SP-8233

SP-8233

SP-8233

Curved blunt tip, 2mm
Angled blunt tip, 4mm
Angled sharp tip, 5mm
0.5mm, 25G
angled sharp tip, 0.6mm, 23G
angled sharp tip, 0.8mm, 20G


Slade LASIK Wire Speculum



CALIPERS

Braunstein Fixed Caliper

Fine tips mark sclerotomy sites for vitreous surgery



Castroviejo Caliper

0-20mm measurements in 0.5mm increments. Scale reading on one side



Luntz-Dodick Trabeculectomy Punch



SP-8233

SP-8233

SP-8233

SP-8233

0.75 mm
1.00 mm
1.20 mm
1.50 mm

LASIK

Crafted with Care

Katzen LASIK Flap Unzipper



SP-8233

double-ended (modified Sinskey hook / bi-convex spatula)

LASIK Multi-purpose Forceps

Combines the features of a spatula, dissector and forceps




SP-8233

Curved jaws, 115mm long

Buratto LASIK Cannula

One front and two side ports. 6 ports for multi-directional irrigation



SP-8233

25G

Serrated Forceps, VR



SP-8233

1.25mm jaws, 23G



Straight toothed Forceps, VR



SP-8233

1.25mm jaws, 23G, VR



Gripping Forceps, VR



SP-8233

1.25mm jaws, 23G, VR



Angled Forceps, membrane peeler, VR



SP-8233

1.25mm, 23G, VR



Angled Forceps, membrane peeler, VR



SP-8233

1.25mm, 23G, VR



Fine Gripping Forceps, VR



SP-8233

3mm jaws, VR, 23G



Side Gripping Forceps



SP-8233

1.25 jaws, 23G



Foreign Body Forceps



SP-8233

Cap jaws, 23G



Micro Capsulorhexis Forceps



SP-8233

1.25 jaws, Curved, 23G



Gripping Forceps, straight, VR



SP-8233

1.25mm jaws, 23G



Gripping Forceps, straight, VR



SP-8233

1.25mm jaws, 23G



End Gripping Forceps



SP-8233

1.25mm jaws, 23G



End Gripping Forceps



SP-8233

3mm jaws, 23G



Flat Gripping Forceps, VR with Tungsten Carbide coating



SP-8233

1.25mm jaws, 23G



Pick Forceps, VR



SP-8233

3mm jaws, 23G



Pointed Gripping Forceps, VR



SP-8233

1.5mm jaws, pointed end, 23G



Offset Gripping Forceps, VR



SP-8233

0.9mm jaws, 23G



ILM Gripping Forceps, VR



SP-8233



2.2mm jaws, 23G, VR

Fine Gripping Forceps, VR



SP-8233



1.25mm jaws, 23G

Serrated Forceps



SP-8233

SP-8233

SP-8233



1.5mm jaws, 23G, VR
1.5mm jaws, curved, 23G, VR
3.0mm jaws, 23G, VR

Gripping Forceps, VR



SP-8233



2.75mm jaws, 23G

Straight Scissors, VR



SP-8233



1.25mm blade, 23G

Straight Scissors, VR



SP-8233



Blunt-tipped, 1.25mm blade, 23G

Curved Scissors, VR



SP-8233



Sharp-tipped, 1.25mm blade, 23G

Angled Scissors, VR



SP-8233



45° angled, 1.5mm blade, 23G

Vertical Scissors, VR



SP-8233



23G

Vertical Scissors, VR



SP-8233



23G

Straight Scissors, VR



SP-8233

SP-8233



Sharp-tipped, 2.25mm blade, 23G
Blunt-tipped, 2.25mm blade, 23G

Curved Scissors, VR



SP-8233



Sharp-tipped, 2.25mm blades, 23G

Straight Scissors, VR



SP-8233



3mm blade, 23G, VR

Vertical Scissors, VR



SP-8233



45° angled, 0.6mm blade, 23G, VR

Vertical Scissors, VR



SP-8233

SP-8233



Curved Scissors, 2.5mm blade, 15mm curved shaft, 23G, VR
Curved Scissors, 3.0mm blade, 25mm curved shaft, 23G, VR

Angled Scissors, VR



SP-8233



45° angled 2mm blade, 23G, VR

Straight Scissors

SP-8233

(Pictures Non Available)
2.5mm blade, 25mm curved shaft, 23G, VR

Irrigation Tips

SP-8233

straight, 0.3mm side port, 21G

SP-8233

45° angled, 0.3mm side port, 21G

SP-8233

Akahoshi mini curved I/A tip, 0.3mm side port, S-shape, 23G

Aspiration Tips

SP-8233

straight, 0.3mm side port, 21G

SP-8233

45° angled, 0.3mm side port, 21G

SP-8233

Akahoshi mini curved I/A tip, 0.3mm side port, S-shape, 23G

Bi-Manual Irrigation Handpieces



SP-8233

Curved, 21G, 0.3mm dual side ports

SP-8233

Curved, 23G, 0.3mm dual side ports

SP-8233

Curved, 21G, 0.5mm dual side ports

SP-8233

Curved, 21G, 0.5mm front port

SP-8233

Curved, 21G, 30° beveled tip

SP-8233

Straight, 23G, 0.4mm dual side ports

SP-8233

Straight, 21G, 0.5mm dual side ports

Bi-Manual Aspiration Handpieces



SP-8233

Curved, 20G, 0.2mm side-port, with capsule polisher

SP-8233

Curved, 20G, 0.3mm side-port, with capsule polisher

SP-8233

Curved, 20G, 0.3mm side-port, cortex extractor with capsule polisher

SP-8233

Curved, 20G, 0.3mm dual side-port, 1.4mm incision

SP-8233

Straight, 20G, 0.3mm dual side-port, 1.4mm incision

SP-8233

Curved, 21G, 0.3mm side-port

SP-8233

Curved, 21G, 0.4mm dual side-port

SP-8233

Curved, 23G, 0.3mm

SP-8233

Curved, 23G, 0.3mm side-port

SP-8233

Curved, 23G, 0.3mm side-port, with polisher

SP-8233

Curved, 23G, 0.4mm side-port

SP-8233

Curved, 23G, 0.4mm side-port, with polisher

Handles

SP-8233

Irrigating Handle

SP-8233

Aspirating Handle

SP-8233



Infusion Handle, M/M connectors

SP-8233

Infusion Handle, M/F connectors

SP-8233

Infusion Handle, M/tubing connectors

	DESCRIPTION	BASIC	FULL
	Utrata-Style Capsulorhexis Forceps, curved (see Code TT0117E at page 6, for complete specs)	✓	✓
	Barraquer Needle Holder, curved (see Code TT0210B at page 9, for complete specs)	✓	✓
	Vannas Capsulotomy Scissors (see Code TT0303H at page 11, for complete specs)	✓	✓
	McPherson Tying Forceps (see Code TT0102D at page 3, for complete specs)	✓	✓
	Kratz-Barraquer Wire Speculum (see Code TT0500H at page 16, for complete specs)	✓	
	Phaco Chopper (see Code TT0430C at page 14, for complete specs)		✓
	Lieberman Speculum, K-wire (see Code TT0500F at page 16, for complete specs)		✓
	Utrata-Style Capsulorhexis Forceps, curved (see Code TT0306C at page 10, for complete specs)		✓
	Bi-Manual Irrigation Handpiece (see Code T301i at page 25, for complete specs)		✓
	Bi-Manual Aspiration Handpiece (see Code T301a at page 25, for complete specs)		✓
	Sterilization Tray / Medium (see Code T104 at page 27, for complete specs)	✓	
	Sterilization Tray / Large (see Code T105 at page 27, for complete specs)		✓


DESCRIPTION



Hunt Chalazion Forceps
(see Code TT0140H at page 5, for complete specs)



Meyerhoeffer Chalazion Curette
SP-8233 2mm
SP-8233 4mm



Iris Scissors
(see Code TT0301I at page 10, for complete specs)




Colibri Toothed Forceps
(see Code TT0111D at page 4, for complete specs)



Sterilization Tray / Medium
(see Code T104 at page 27, for complete specs)

STERILIZATION TRAYS



Poly-Carbonate Sterilization Trays with Silicone Mat

Code	Sizes / Dimensions	
SP-8233	Small	(16cm x 8cm x 2cm)
SP-8233	Medium	(20cm x 13cm x 2cm)
SP-8233	Large	(24cm x 14cm x 2cm)