


VETERINARY ORTHOPEDICS IMPLANTS



MANUFACTURER OF VETERINARY ORTHOPEDICS IMPLANTS

- 1.5 mm Locking Cuttable Plates
- 2 mm Locking Cuttable Plates
- 2.4 mm Locking Cuttable Plates
- 2.7 mm Locking Cuttable Plates
- 3.5 mm Locking Cuttable Plates

Each Plate has 20 holes

Plate can cut according to required length

available in Titanium



PROFILE

Items	1.5 mm Plate	2 mm Plate	2.4 mm Plate	2.7 mm Plate	3.5 mm Plate
Thickness	1.5	2	2.5	2.5	3
Width	4	5	6	7	9
Hole Spacing	5	6	7	9	11
Screw	1.5 mm Locking	2 mm Locking	2.4 mm Locking	2.7 mm Locking	3.5 mm Locking

CODES

Size (mm)	Holes	Length (mm)	TI (codes)
1.5	20	100	LCPS15
2	20	120	LCPS20
2.4	20	140	LCPS24
2.7	20	180	LCPS27
3.5	20	220	LCPS35



2 mm Locking Cuttable T Plates
2.4 mm Locking Cuttable T Plates
2.7 mm Locking Cuttable T Plates
3.5 mm Locking Cuttable T Plates

Each Plate has 7 holes
Plate can cut according to required length

available in Titanium



PROFILE

Items	2 mm Plate	2.4 mm Plate	2.7 mm Plate	3.5 mm Plate
Thickness	2	2.5	2.5	3
Width	6	7	8	10
Total Length	47	55	62	77
Screw	2 mm Locking	2.4 mm Locking	2.7 mm Locking	3.5 mm Locking


CODES

Size (mm)	Holes	Length (mm)	TI (codes)
2	7	47	LCTP20
2.4	7	55	LCTP24
2.7	7	62	LCTP27
3.5	7	77	LCTP35



Left & Right Side

available in Titanium



PROFILE

Items	2 mm Plate	2.4 mm Plate	2.7 mm Plate	3.5 mm Plate
Thickness	1.5	2	2	2.5
Width	6	7	8	10
Hole Spacing	8	9	11	13
Screw	2 mm Locking 2 mm Cortical	2.4 mm Locking 2.4 mm Cortical	2.7 mm Locking 2.7 mm Cortical	3.5 mm Locking 3.5 mm Cortical 4 mm Cancellous

CODES

	Left Codes				Right Codes			
Holes	2 mm	2.4 mm	2.7 mm	3.5 mm	2 mm	2.4 mm	2.7 mm	3.5 mm
4	LSP01L	LSP06L	LSP011L	LSP016L	LSP01R	LSP06R	LSP011R	LSP016R
5	LSP02L	LSP07L	LSP012L	LSP017L	LSP02R	LSP07R	LSP012R	LSP017R
6	LSP03L	LSP08L	LSP013L	LSP018L	LSP03R	LSP08R	LSP013R	LSP018R
7	LSP04L	LSP09L	LSP014L	LSP019L	LSP04R	LSP09R	LSP014R	LSP019R
8	LSP05L	LSP10L	LSP015L	LSP020L	LSP05R	LSP10R	LSP015R	LSP020R



2.4 mm Locking TPLO Plate

2.7 mm Locking TPLO Plate

3.5 mm Locking TPLO Plate

Left & Right Side

available in Titanium



PROFILE

Items	2.4 mm Plate	2.7 mm Plate	3.5 mm Plate	3.5 mm Broad Plate
Thickness	2.5	2.5	3.5	3.5
Width	7	8	10	13
Total Length	50/58	55/64	67/77	75
Screw	2.4 mm Locking 2.4 mm Cortical	2.7 mm Locking 2.7 mm Cortical	3.5 mm Locking 3.5 mm Cortical 4 mm Cancellous	3.5 mm Locking 3.5 mm Cortical 4 mm Cancellous

CODES

Holes	Left Codes			Right Codes		
	2.4 mm	2.7 mm	3.5 mm	2.4 mm	2.7 mm	3.5 mm
3	TPLO1L	TPLO3L	TPLO5L	TPLO1R	TPLO3R	TPLO5R
4	TPLO2L	TPLO4L	TPLO6L	TPLO2R	TPLO4R	TPLO6R



2 mm Locking Compression Plates
2.4 mm Locking Compression Plates
2.7 mm Locking Compression Plates
3.5 mm Locking Compression Plates
5 mm Locking Narrow Compression Plates
5 mm Locking Broad Compression Plates

available in Titanium & SS 316L



HOLES

4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18

PROFILE

Items	2 mm Plate	2.4 mm Plate	2.7 mm Plate	3.5 mm Plate	5 mm Narrow Plate	5 mm Broad Plate
Thickness	2	2.5	2.5	3.5	4.5	6
Width	6	7	9	11	13.5	17.5
Hole Spacing	8	9	11	13	18	18
Screw	2 mm Locking 2 mm Cortical	2.4 mm Locking 2.4 mm Cortical	2.7 mm Locking 2.7 mm Cortical	3.5 mm Locking 3.5 mm Cortical 4 mm Cancellous	5 mm Locking 4.5 mm Cortical 5 mm Cancellous	5 mm Locking 4.5 mm Cortical 5 mm Cancellous

2 mm Dynamic Compression Plates
2.4 mm Dynamic Compression Plates
2.7 mm Dynamic Compression Plates

3.5 mm Dynamic Compression Plates
4.5 mm Narrow Dynamic Compression Plates
4.5 mm Broad Dynamic Compression Plates

available in Titanium & SS 316L



HOLES

4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18

PROFILE

Items	2 mm Plate	2.4 mm Plate	2.7 mm Plate	3.5 mm Plate	4.5 mm Narrow Plate	4.5 mm Broad Plate
Thickness	2	2.5	2.5	3.5	4	5
Width	6	7	8	10	12	16
Hole Spacing	6	7	8	12	16	16
Centre Hole Spacing	8	10	12	16	25	25
Screw	2 mm cortical	2.4 mm cortical	2.7 mm cortical	3.5 mm cortical	4.5 mm cortical	4.5 mm cortical



Left & Right Side

available in Titanium



Size	Thickness (mm)	Direction	Length (mm)	TI (codes)
2/2.7 mm	2	Left	90	LPAPA1L
2/2.7 mm	2	Right	90	LPAPA1R
2.7/3.5 mm	2.5	Left	120	LPAPA2L
2.7/3.5 mm	2.5	Right	120	LPAPA2R

2/2.4 mm Locking Pancarpal Arthrodesis Plates
2.4/2.7 mm Locking Pancarpal Arthrodesis Plates
2.7/3.5 mm Locking Pancarpal Arthrodesis Plates

available in Titanium & SS 316L



Size (mm)	Thickness (mm)	Holes	Length (mm)	Codes
2/2.4	2.5	9	90	LPAP1
2/2.4	2.5	11	97	LPAP2
2.4/2.7	2.5	9	97	LPAP3
2.4/2.7	2.5	11	115	LPAP4
2.7/3.5	3	9	110	LPAP5
2.7/3.5	3	11	120	LPAP6

2 mm Locking T Plates Clover Leaf
2.4 mm Locking T Plates Clover Leaf

2.7 mm Locking T Plates Clover Leaf
3.5 mm Locking T Plates Clover Leaf

available in Titanium



PROFILE

Items	2 mm Plate	2.4 mm Plate	2.7 mm Plate	3.5 mm Plate
Thickness	1.5	2	2	2.5
Width	6	7	8	10
Hole Spacing	8	9	11	13
Screw	2 mm Locking 2 mm Cortical	2.4 mm Locking 2.4 mm Cortical	2.7 mm Locking 2.7 mm Cortical	3.5 mm Locking 3.5 mm Cortical 4 mm Cancellous

CODES

Holes	2 mm	2.4 mm	2.7 mm	3.5 mm
3	LTPCL01	LTPCL06	LTPCL11	LTPCL16
4	LTPCL02	LTPCL07	LTPCL12	LTPCL17
5	LTPCL03	LTPCL08	LTPCL13	LTPCL18
6	LTPCL04	LTPCL09	LTPCL14	LTPCL19
8	LTPCL05	LTPCL10	LTPCL15	LTPCL20



- 1.5 mm Locking Double T Plate
- 2.0 mm Locking Double T Plate
- 2.4 mm Locking Double T Plate
- 2.7 mm Locking Double T Plate
- 3.5 mm Locking Double T Plate

Each Plate has 20 holes
One side three head holes and other side two head holes.
Plate can cut according to required length

available in Titanium



PROFILE

Items	1.5 mm Plate	2 mm Plate	2.4 mm Plate	2.7 mm Plate	3.5 mm Plate
Thickness	1.5	2	2.5	2.5	3
Width	4	6	6.5	7	9
Hole Spacing	4	6	6.5	7	9
Screw	1.5 mm Locking	2 mm Locking	2.4 mm Locking	2.7 mm Locking	3.5 mm Locking

CODES

Size (mm)	Holes	Length (mm)	TI (codes)
1.5	20	130	LDTP15
2	20	155	LDTP20
2.4	20	170	LDTP24
2.7	20	185	LDTP27
3.5	20	235	LDTP35



Any length can be manufacture

available in Titanium & SS 316L

2 mm Cortical Screws



2.4 mm Cortical Screws



2.7 mm Cortical Screws



3.5 mm Cortical Screws



4.5 mm Cortical Screws



6/ 7/ 8 mm Veterinary Nails Straight



Nail Length 250 mm


Nail can cut according to required length

3.4 mm Screws are used

available in Titanium & SS 316L

3.4 mm Nail Screws





7/ 8/9 mm Veterinary Nails Angled

Nail Length

140 / 160 / 180 / 200 / 220 / 240/ mm

3.4 mm Screws are used

available in Titanium & SS 316L

3.4 mm Nail Screws



Toggle Pins

available in Titanium



Size (mm)	Hole Dia (mm)	Length (mm)	TI (codes)
3	1.2	10	TGP1
3.5	1.5	15	TGP2
4	2	20	TGP3

Suture Buttons (ACL Endo Button)

available in Titanium



Thickness : 2 mm
 Width : 6 mm
 length : 10 mm
 Hole dia : 2.5 mm



Thickness : 2 mm
 Width : 4 mm
 length : 10 mm
 Hole dia : 1.5 mm



Thickness : 2 mm
 Width : 5 mm
 length : 16 mm
 First & last hole dia : 1.5 mm
 Two middle hole dia : 2 mm



PROFILE

Size	Connecting Rod Dia (mm)	Pin Size (mm)
Mini	3	1.5 1.8
Small	3.5	2.5 3.0
Large	3.5	3.2 3.5



DIA 3 mm | 3.5 mm

LENGTH 100 mm | 125 mm | 150 mm | 200 mm | 250 mm | 300 mm



DIA 1.5 mm | 1.8 mm | 2.5 mm | 3.0 mm | 3.2 mm | 3.5 mm

LENGTH 80 mm | 100 mm | 125 mm | 150 mm | 200 mm | 300 mm

3.5 / 4 / 5 / mm Locking Cancellous Screws


Any length can be manufacture

available in Titanium & SS 316L

3.5 mm Locking Cancellous Screws




4 mm Locking Cancellous Screws



5 mm Locking Cancellous Screws, Part Thread



5 mm Locking Cancellous Screws, Full Thread



2.7 / 3.5 / 5 / mm Cannulated Screws


Any length can be manufacture

available in Titanium & SS 316L

2.7 mm Cannulated Screws



3.5 mm Cannulated Screws



5 mm Cannulated Screws



2 / 2.7 / 3.5 / 4.5 / 5.5 mm Spine Screws

available in Titanium

Length (mm): 14 mm - 60 mm

MonoAxial Screws




MonoAxial Reduction Screws



PolyAxial Screws



PolyAxial Reduction Screws



Crosslink Connector / Transverse Connector



Length (mm) : 80 | 90 | 100 | 120 | 150 | 200

Crosslink Connector / Transverse Connector



Length (mm) : 80 | 100 | 120 | 140 | 160 | 180 | 200 | 250 | 300

Single
Trocár



Double
Trocár



One Side
Thread



Fully
Thread



Sizes

Dia (mm) : 0.8 to 6

Length (mm) : 80 to 300

Crafted with Care

V159 Drill Bits

Cannulated Drill Bit



Cannulated Drill Bit, Quick Coupling



Drill Bit



Drill Bit, Quick Coupling



Sizes

Dia (mm) : 1 mm - 5 mm

Length (mm) : 60 to 250



Sizes

16 SWG (1.60 mm)

18 SWG (1.25 mm)

20 SWG (0.90 mm)

22 SWG (0.70 mm)

24 SWG (0.55 mm)


26 SWG (0.45 mm)

28 SWG (0.35 mm)

30 SWG (0.30 mm)



V001 Screw Driver - Torque



Available for
2.7 / 3.5 / 5 Locking Screws

V002 Screw Driver - Fibre Handle

Crafted with Care



Available for
2 / 2.4 / 2.7 / 3.5 / 5 Locking Screws
1.5 / 2 / 2.4 / 2.7 / 3.5 / 4.5 Cortical Screws
2.7 / 3.5 / 5 / Cancellous Screws

V003 Screw Driver - Holding Sleeve



Available for
2.4 / 2.7 / 3.5 / 5 Locking Screws
2 / 2.4 / 2.7 / 3.5 / 4.5 Cortical Screws
2.7 / 3.5 / 5 / Cancellous Screws

V004 Screw Driver - With Sleeve



Available for
1.5 / 2 Locking Screws
1.5 / 2 Cortical Screws

Crafted with Care

SURGIPORTAL



Available for
2.7 / 3.5 / 5 Cannulated Screws

V010 Bone Tap T Handle



Available for

Cortical Screws: 1.5 | 2 | 2.4 | 2.7 | 3.5 | 4.5

Cancellous Screws: 3.5 | 5

Cannulated Screws: 2.7 | 3.5 | 5

V011 Bone Tap Coupling Shaft



Available for

Cortical Screws: 1.5 | 2 | 2.4 | 2.7 | 3.5 | 4.5

Cancellous Screws: 3.5 | 4 | 4.5 | 6.5 | 7 | 7.3

Cannulated Screws: 2.5 | 3.5 | 4 | 4.5 | 6.5 | 7 | 7.3

V012 Countersink Coupling Shaft



Available for

Cortical Screws: 2.4 | 2.7 | 3.5 | 4.5

Cancellous Screws: 2.7 | 3.5 | 4.5

Cannulated Screws: 2.7 | 3.5 | 4.5

V013 Countersink T Handle



Available for

Cortical Screws: 2.4 | 2.7 | 3.5 | 4.5

Cancellous Screws: 3.5 | 4.5

Cannulated Screws: 2.7 | 3.5 | 4.5



Available for

2 / 2.4 / 2.7 / 3.5 / 5 Locking Screws

1.5 / 2 / 2.4 / 2.7 / 3.5 / 4.5 Cortical Screws

3.5 / 4 / 5 / 6.5 Cancellous Screws



V008 Screw Driver Coupling Handle, Plastic



V009 Screw Driver - PolyAxial Screw



Available for : 2 / 2.4 / 2.7 / 3.5 / 4.5 / 5.5 / 6.5 PolyAxial Screws

V018 Self Centering Bone Holding Forceps Richat Lock

Sizes (inches)
6 | 7 | 8 | 9 | 10 | 11 | 12



V019 Self Centering Bone Holding Forceps
Speed Lock



V020 Hey Groves Bone Holding Forceps

Sizes (inches)
6 | 8 | 10 | 12




Length (inches)
7 | 8 | 9.5 | 10.5



V022 Patella Bone Holding Forceps 4 Prong

Length (inches)
9.5 | 10.5



V023 Lowman Bone Clamp





Length (inches) : 6




V027 Bone and Plate Holding Forceps - Type 1



Length (inches) :






Length (inches) : 9.5 | 11







Size: 7.5 inches



Size: 7.5 inches



V042 Kerrison Bone Punch



V043 Meniscal knife



V044 Osteotome Straight/Curved Fibre Handle Length : 10 inches



Tip Sizes (mm) for Straight
5 | 10 | 15 | 20 | 25 | 30

Tip Sizes (mm) for Curved
5 | 10 | 15 | 20 | 25 | 30


V045 Osteotome Straight/Curved SS Handle Length : 8 inches



Tip Sizes (mm) for Straight
5 | 10 | 15 | 20 | 25 | 30

Tip Sizes (mm) for Curved
5 | 10 | 15 | 20 | 25 | 30

V046 Bone Gouge Fiber Handle Length : 10 inches



Tip Sizes (mm) for Straight
5 | 10 | 15 | 20 | 25 | 30

Tip Sizes (mm) for Curved
5 | 10 | 15 | 20 | 25 | 30

V047 Bone Gouge SS Handle Length: 8 inches



Tip Sizes (mm) for Straight
5 | 10 | 15 | 20 | 25 | 30

Tip Sizes (mm) for Curved
5 | 10 | 15 | 20 | 25 | 30

V048 Capener Bone Gouge Size: 10 inches
Sizes (mm) for Straight
5 | 10 | 15 | 20 | 25 | 30



V049 Chisel Bone Gouge Size: 10 inches



Tip Sizes (mm) for Straight
5 | 10 | 15 | 20 | 25 | 30

Tip Sizes (mm) for Curved
5 | 10 | 15 | 20 | 25 | 30

V052 Gelpi Retractor

Length (inches)
6 | 8 | 10




Length (inches)
7 | 9



V054 Lamina Spreaders

Length (inches)
7 | 9





Crafted with Care
SURGIPORTAL




Length : 8 inches

Tip Sizes (mm) : 5 | 10 | 15 | 20 | 25 | 30



Length : 8 inches

Tip Sizes (mm) : 5 | 10 | 15 | 20 | 25 | 30



Length : 10 inches

Tip Sizes

Small | Medium | Large



V069 Bristow Elevator

Length : 10 inches





Length : 9 inches

Jaw Sizes : Small | Medium | Large



V077 Bone Curette Double Ended



Length : 9 inches

Cup Sizes : Small | Medium | Large

V078 Bone Curette Fiber Handle



Cup Sizes (mm)

3 | 4 | 5 | 6 | 9 | 12

Length : 9 inches

V079 Bone Curette SS Handle



Cup Sizes (mm)

3 | 4 | 5 | 6 | 9 | 12

Length : 15 inches

V080 Ring Curette 15 inches



Length : 15 inches

V073 Bone File Flat



Length : 8 inches

End Sizes : Small | Medium | Large

V074 Bone File Flat Extended



Length (inches) : 6 | 8 | 10

V075 Bone File Round



Length : 8 inches

End Sizes : Small | Medium | Large

V076 Bone File Round Extended



End Sizes : Small | Medium | Large

V081 Bone Hook



Length : 10 inches

End Sizes : Small | Medium | Large

V082 Bone Hook T Handle



Length : 8 inches

End Sizes : Small | Medium | Large



Length : 6 inches


V084 Mensical Probe



V085 Hammer SS Face

Weight (Grams) : 250 | 350 | 500

Length : 10 inches



V086 Hammer Nylon and SS Face

Weight (Grams) : 250 | 350 | 500


Length : 10 inches



V087 Hammer Nylon Face both Side

Weight (Grams) : 250 | 350 | 500

Length : 10 inches






Crafted with Care

SURGIPORTAL





V097 Gigli saw handle set of 2



K Wire




Sizes (Inches)
6 | 8 | 10 | 12 | 14 | 16 | 18 | 20

V099 Wire Passer Fibre Handle



Sizes
Small | Medium | Large

V100 Wire Passer SS Size : 11 inches





Length : 4 inches



Length : 10 inches



V105 Rod Bender

Length : 15 inches



V106 Plates Bender Roller

Length : 12 inches




V107 Plate Bending Press Bench



V108 Plate Bending Template

Available Sizes
Small | Medium | Large





Size (mm): 1.5/2.0 | 2.0/2.7 | 2.5/3.5 | 3.2/4.5 | 3.2/6.5

V129 Drill Guide Neutral and Load



Available Sizes: 1.5/2.0 | 2.0/2.7 | 2.5/3.5 | 3.2/4.5

1.5 MM Drill for 2.0 mm Screws | 2.0 MM Drill for 2.7 mm Screws

2.5 MM Drill for 3.5 mm Screws | 3.2 MM Drill for 4.5 mm Screws



Available for : 2 / 2.4 / 2.7 / 3.5 / 5 mm Locking Screws

V131 Drill Bit Sleeve



Available for following Drills (in mm) : 2 | 2.5 | 3 | 3.5 | 4 | 4.5

V132 Tap Sleeve



Available for following Taps (in mm)

1.5 | 2 | 2.5 | 3 | 3.5 | 4 | 4.5

V126 BP Handle

Sizes: 3 | 4



V127 Tissues & Disc Forcep

Available Sizes for Straight (in mm)
2 | 3 | 4



Available Sizes for Upward (in mm)
2 | 3 | 4



Available Sizes for Downward (in mm)
2 | 3 | 4





30 mm Measuring




60 mm Measuring



110 mm Measuring



Available Sizes
Small | Medium | Large





V141 Drill Chuck & Key

Available Sizes
Small | Medium | Large




V142 Schanz Pin Introducer

Length : 4 inches



V143 Steinmann Pin Holder

Length (inches)
6 | 8 | 10 | 12 | 14



V144 Pin Punch



Pin Punch, 6 mm dia, 4 inches length
Pin Punch, 6 mm dia, 5.5 inches length



Length : 6 inches

Available for following Screws
2.4 | 2.7 | 3.5 | 4.5 | 5



V148 Broken Screw Head Removal Ratchet Lock



Length : 10 inches

V149 Broken Screw Head Removal Shear

Length : 8 inches



Available for following Screws

2.7 | 3.5 | 4.5

Length : 6 inches



V151 Screw Holding Forceps

Length : 3 inches



V152 Bandage Cutting Scissor

Length : 6 inches



V153 Esmarch Bandage


Sizes (inches)

3 | 4 | 6



V154 Screws Washer Hanger





Battery Operated Oscillation Saw with


1 Machine | 2 Batteries | 3 Blades
1 Charger | 1 Sterilization Ring | 1 Carry Case






Drill Machine Battery Operated (cordless)

1 Machine | 2 Batteries | 1 Charger
1 Key | 1 Chuck | 1 Carry Case





2.4 mm Threaded Drill Guide for 1.8 mm Drill Bit

2.7 mm Threaded Drill Guide for 2 mm Drill Bit



Depth Gauge for 2 / 2.4 mm Screws




Depth Gauge for 2.7 mm Screws




Double Drill Sleeve (Universal Drill Sleeve) 1.8/2.4 mm

Double Drill Sleeve (Universal Drill Sleeve) 2.0/2.7 mm



Drill Bit 1.8 x 100 mm, Quick Coupling
Drill Bit 2 x 100 mm, Quick Coupling
Drill Bit 2.4 x 100 mm, Quick Coupling
Drill Bit 2.7 x 100 mm, Quick Coupling






Drill Bit 2.5 mm x 125 mm, Quick Coupling
Drill Bit 3.5 mm x 130 mm, Quick Coupling







Drill Guide 3.5 mm for Neutral and Load Position



Periosteal Elevator Fiber Handle 10 mm



Reduction Forcep Serrated Speed Lock 190 mm



Screw Extraction T Handle for 3.5 mm Screws





Tap T Handle for 3.5 mm Cortical Screws
Tap T Handle for 4 mm Cancellous Screws








Crafted with Care




Trephine T Handle









Drill Bit 3.2 mm x 130 mm Quick Coupling

Drill Bit 4.5 mm x 130 mm Quick Coupling







Periosteal Elevator Fiber Handle 5 mm



Periosteal Elevator Fiber Handle 10 mm










Crafted with Care

SURGIPORTAL









Screw Driver for Offset Screw (Nut)



Screw Holder for Offset Screw (Nut)



PolyAxial Screw Driver 2 Piece



Screw Holder for Offset Screw (Nut)






Spine Rod Pusher Curved




Spine Rod Pusher Straight









Crafted with Care


S.No.	Description	
1	Jig for Nail	1
2	Conical Bolt for Jig	1
3	Protection Sleeve	1
4	Trocar	1
5	Drill Sleeve	1
6	Depth Gauge measuring upto 110 mm	1
7	Screw Driver	1
8	Bone Awl Cannulated	1
9	Spanner	1
10	Rod Handle	1
11	Rod	1
12	Ram	1
13	2 mm Guide Wire	2
14	Drill Bit, 2.7 mm x 225 mm	2
15	Reamer T handle Cannulated 6 mm	1
16	Reamer T handle Cannulated 7 mm	1
17	Reamer T handle Cannulated 8 mm	1
18	Graphics Aluminum Box with Silicone Fittings	1



Depth Gauge measuring upto 110 mm








Crafted with Care

S.No.	Description	
1	Jig for Nail	1
2	Conical Bolt for Jig	1
3	Protection Sleeve	1
4	Trocar	1
5	Drill Sleeve	1
6	Depth Gauge measuring upto 110 mm	1
7	Screw Driver	1
8	Bone Awl Cannulated	1
9	Spanner	1
10	Rod Handle	1
11	Rod	1
12	Ram	1
13	2 mm Guide Wire	2
14	Drill Bit, 2.7 mm x 225 mm	2
15	Reamer T handle Cannulated 6 mm	1
16	Reamer T handle Cannulated 7 mm	1
17	Reamer T handle Cannulated 8 mm	1
18	Graphics Aluminum Box with Silicone Fittings	1

Drill Bit, 2.7 mm x 225 mm



Reamer T handle Cannulated 6 mm




Crafted with Care


SURGIPORTAL



SURGI PORTAL








Depth Gauge measuring upto 110 mm





S.NO.	DESCRIPTION	QTY
1	3.5 mm Cortical Screws	140
2	3.5 mm Locking Screws	140
3	3.5/4 mm Cancellous Screws	92



S.NO.	DESCRIPTION	QTY
1	2 mm Locking Screws	48
2	2 mm Cortical Screws	48
3	2.4 mm Locking Screws	48
4	2.4 mm Cortical Screws	48
5	2.7 mm Locking Screws	48
6	2.7 mm Cortical Screws	48



S.NO.	DESCRIPTION	QTY
1	5 mm Locking Screws	130
2	4.5 mm Cortical Screws	104
3	6.5 mm Cancellous Screws	98







S.NO.	DESCRIPTION	QTY
1	2.7 mm Cannulated Screws	33
2	3.5 mm Cannulated Screws	54
3	4.0 mm Cannulated Screws	54



S.NO.	DESCRIPTION	QTY
1	2.0 mm Polyaxial Screws	28
2	2.7 mm Polyaxial Screws	28
3	3.5 mm Polyaxial Screws	56
4	4.5 mm Polyaxial Screws	28
5	5.5 mm Polyaxial Screws	28



Company Name: SurgiPortal

Phone: +92 339 0059380

Email: info@surgiportal.com

Location: Ikram Plaza, Near KFC, Kashmir Road, Sialkot, Pakistan